

Anti-LTA antibody [55] (STJ16100772)

STJ16100772

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-LTA is suitable for use in Immunohistochemistry, Flow Cytometry, ELISA, Immunofluorescence and Western Blot research applications.
Applications	IHC-F/FC/ELISA/IF/WB
Host/Source	Mouse
Reactivity	

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	55
Concentration	> 200 ug/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.1% BSA and 0.02% sodium azide
Isotype	IgG3
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID	4049
Gene Symbol	LTA
Uniprot ID	TNFB_HUMAN
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	