

Anti-CD46 antibody [M177] (STJ16100774)

STJ16100774

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD46 is suitable for use in Immunohistochemistry, Flow Cytometry, ELISA, Immunoprecipitation and Western Blot research applications.
Applications	IHC-F/FC/ELISA/IP/WB
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	M177
Concentration	100 ug/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.1% BSA and 0.02% sodium azide
Isotype	IgG1
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID	4179
Gene Symbol	CD46
Uniprot ID	MCP_HUMAN
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	