

**Anti-PR3 antibody [PRG-2] (STJ16100777)**

STJ16100777

**GENERAL INFORMATION**

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Mouse monoclonal antibody anti-PR3 is suitable for use in Flow Cytometry, ELISA and Western Blot research applications.
<b>Applications</b>	FC/ELISA/WB
<b>Host/Source</b>	Mouse
<b>Reactivity</b>	Human

**PRODUCT PROPERTIES**

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	PRG-2
<b>Concentration</b>	100 ug/mL
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	
<b>Dilution Range</b>	
<b>Formulation</b>	PBS with 0.1% BSA and 0.02% sodium azide
<b>Isotype</b>	IgG1
<b>Storage Instruction</b>	Store at 2-8°C upon receipt.

**TARGET INFORMATION**

<b>Gene ID</b>	<a href="#">5657</a>
<b>Gene Symbol</b>	<a href="#">PRTN3</a>
<b>Uniprot ID</b>	<a href="#">PRTN3_HUMAN</a>
<b>Immunogen</b>	
<b>Immunogen Region</b>	
<b>Specificity</b>	
<b>Immunogen Sequence</b>	