

Anti-Lactoferrin antibody [265-1K1] (STJ16100779)

STJ16100779

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Lactoferrin is suitable for use in Immunohistochemistry, ELISA and Western Blot research applications.
Applications	IHC-F/ELISA/WB
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	265-1K1
Concentration	100 ug/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.1% BSA and 0.02% sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	4057
Gene Symbol	LTF
Uniprot ID	TRFL_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen	
Sequence	