

Anti-ASGPR antibody [8D7] (STJ16100783)

STJ16100783

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-ASGPR is suitable for use in Immunohistochemistry, Flow Cytometry, ELISA, Immunofluorescence and Western Blot research applications.
Applications	IHC-F/FC/ELISA/IF/WB
Host/Source	Mouse
Reactivity	Rat/Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	8D7
Concentration	100 ug/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.1% BSA and 0.02% sodium azide
Isotype	IgG1
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen
Region
Specificity
Immunogen
Sequence