

Anti-HSV2 antibody [T303] (STJ16100789)

STJ16100789

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Herpes simplex virus type 2 is suitable for use in Western Blot, Immunocytochemistry and ELISA research applications.
Applications	WB/ICC/ELISA
Host/Source	Mouse
Reactivity	Herpes Simplex Virus

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	T303
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Purified by protein-A affinity chromatography.
Dilution Range	
Formulation	Phosphate buffered saline (PBS) solution with 15 mM sodium azide
Isotype	IgG2a
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Extract of HSV-2 infected VERO (green monkey kidney) cells.
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	