

**Anti-Tenascin-C antibody [EB2] (STJ16100797)**

STJ16100797

**GENERAL INFORMATION**

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Mouse monoclonal antibody anti-Tenascin-C is suitable for use in Immunohistochemistry, Western Blot and ELISA research applications.
<b>Applications</b>	IHC-F/WB/ELISA
<b>Host/Source</b>	Mouse
<b>Reactivity</b>	Human/Rabbit/Guinea Pig

**PRODUCT PROPERTIES**

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	EB2
<b>Concentration</b>	100 ug/mL
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	
<b>Dilution Range</b>	
<b>Formulation</b>	PBS with 1% BSA & 0.1% sodium azide
<b>Isotype</b>	IgG1
<b>Storage Instruction</b>	Store at 2-8°C upon receipt.

**TARGET INFORMATION**

<b>Gene ID</b>	3371
<b>Gene Symbol</b>	TNC
<b>Uniprot ID</b>	TENA_HUMAN
<b>Immunogen</b>	Purified tenascin-C from human fetal fibroblasts
<b>Immunogen Region</b>	
<b>Specificity</b>	
<b>Immunogen</b>	
<b>Sequence</b>	