

**Anti-Beta1-Integrin antibody [DF5] (STJ16100803)**

STJ16100803

**GENERAL INFORMATION**

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Mouse monoclonal antibody anti-Beta1-Integrin is suitable for use in Immunohistochemistry, Western Blot and ELISA research applications.
<b>Applications</b>	IHC-P/IHC-F/WB/ELISA
<b>Host/Source</b>	Mouse
<b>Reactivity</b>	Human

**PRODUCT PROPERTIES**

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	DF5
<b>Concentration</b>	100 ug/mL
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	
<b>Dilution Range</b>	IHC-P IHC-F WB 1:1000 ELISA
<b>Formulation</b>	PBS with 1% BSA & 0.1% sodium azide
<b>Isotype</b>	IgG1
<b>Storage Instruction</b>	Store at 2-8°C upon receipt.

**TARGET INFORMATION**

<b>Gene ID</b>	
<b>Gene Symbol</b>	
<b>Uniprot ID</b>	
<b>Immunogen</b>	Wheat germ agglutinin binding glycoproteins of human fibroblasts
<b>Immunogen Region</b>	
<b>Specificity</b>	
<b>Immunogen</b>	
<b>Sequence</b>	