

Anti-Epidermal Growth Factor antibody [E5] (STJ16100805)

STJ16100805

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Epidermal Growth Factor is suitable for use in Immunohistochemistry, Western Blot and ELISA research applications.
Applications	IHC-P/IHC-F/WB/ELISA
Host/Source	Mouse
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	E5
Concentration	5 ug/mL
Conjugation	Unconjugated
Purification	
Dilution Range	IHC-P 1:5 IHC-F 1:5 WB ELISA
Formulation	Culture medium with 0.01M Sodium Azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	EGF coupled to bovine serum albumin and to tetanus toxoid
Immunogen Region	
Specificity	
Immunogen Sequence	