

**Anti-Epidermal Growth Factor antibody [F5] (STJ16100806)**

STJ16100806

**GENERAL INFORMATION**

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Mouse monoclonal antibody anti-Epidermal Growth Factor is suitable for use in Immunohistochemistry, Western Blot and ELISA research applications.
<b>Applications</b>	IHC-P/IHC-F/WB/ELISA
<b>Host/Source</b>	Mouse
<b>Reactivity</b>	Human/Mouse

**PRODUCT PROPERTIES**

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	F5
<b>Concentration</b>	10 ug/mL
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	
<b>Dilution Range</b>	IHC-P 1:20 IHC-F 1:20 WB ELISA
<b>Formulation</b>	Culture medium with 0.01M Sodium Azide
<b>Isotype</b>	IgM
<b>Storage Instruction</b>	Store at 2-8°C upon receipt.

**TARGET INFORMATION**

<b>Gene ID</b>	
<b>Gene Symbol</b>	
<b>Uniprot ID</b>	
<b>Immunogen</b>	mouse EGF coupled to bovine serum albumin and tetanus toxoid
<b>Immunogen Region</b>	
<b>Specificity</b>	
<b>Immunogen Sequence</b>	