

Anti-Bromodeoxyuridine antibody [IIB5] (STJ16100807)

STJ16100807

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Bromodeoxyuridine is suitable for use in Flow Cytometry, Immunohistochemistry, Immunohistochemistry and Immunocytochemistry research applications.
Applications	FC/IHC/IHC-F/IHC-P/ICC
Host/Source	Mouse
Reactivity	Other

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	IIB5
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	
Dilution Range	FC 1:100-1:200
	IHC
	IHC-F
	IHC-P
Formulation	PBS with 0.09% sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	BrdU conjugated to Bovine Serum albumin
Immunogen Region	
Specificity	
Immunogen Sequence	

