

Anti-Chlamydia Trachomatis antibody [3] (STJ16100808)

STJ16100808

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Chlamydia Trachomatis is suitable for use in Immunohistochemistry and Immunofluorescence research applications.
Applications	IHC-P/IF
Host/Source	Mouse
Reactivity	Bacteria

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	3
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	Supernatant with 0.1% sodium azide
Isotype	IgG3
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen
Region
Specificity
Immunogen
Sequence