

Anti-Alpha-IIb-Integrin antibody [CA3] (STJ16100813)

STJ16100813

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Alpha-IIb-Integrin is suitable for use in Immunohistochemistry, Western Blot and ELISA research applications.
Applications	IHC-F/WB/ELISA
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	CA3
Concentration	100 ug/mL
Conjugation	Unconjugated
Purification	
Dilution Range	IHC-F WB 1:1000 ELISA
Formulation	PBS with 1% BSA & 0.1% sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Partially purified platelet glycoprotein IIb/IIIa complex
Immunogen Region	
Specificity	
Immunogen Sequence	