

Anti-Alpha-Fodrin antibody [AA6] (STJ16100815)

STJ16100815

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Alpha-Fodrin is suitable for use in Immunohistochemistry, Western Blot and ELISA research applications.
Applications	IHC-F/WB/ELISA
Host/Source	Mouse
Reactivity	Human/Rat/Chicken/Rabbit/Guinea Pig

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	AA6
Concentration	100 ug/mL
Conjugation	Unconjugated
Purification	
Dilution Range	IHC-F WB 1:1000 ELISA
Formulation	PBS with 1% BSA & 0.1% sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Chicken red blood cell membranes purified by hypotonic lysis and mechanical enucleation
Immunogen Region	
Specificity	
Immunogen Sequence	