

Anti-Cytokeratin 7 antibody [4F5] (STJ16100818)

STJ16100818

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Cytokeratin 7 is suitable for use in Immunohistochemistry and Western Blot research applications.
Applications	IHC-F/WB
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	4F5
Concentration	100 ug/mL
Conjugation	Unconjugated
Purification	
Dilution Range	IHC-F WB 1:1000
Formulation	PBS with 1% BSA & 0.1% sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	3855
Gene Symbol	KRT7
Uniprot ID	K2C7_HUMAN
Immunogen	Cytokeratins isolated from cultured pig kidney epithelial cells
Immunogen Region	
Specificity	
Immunogen Sequence	