

Anti-Erythroid Alpha-Spectrin antibody [AF10] (STJ16100820)

STJ16100820

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Erythroid Alpha-Spectrin is suitable for use in Immunohistochemistry, Western Blot and Immunocytochemistry research applications.
Applications	IHC-F/WB/ICC
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	AF10
Concentration	100 ug/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 1% BSA & 0.1% sodium azide
Isotype	IgG1
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID	6708
Gene Symbol	SPTA1
Uniprot ID	SPTA1_HUMAN
Immunogen	Ghost proteins of human red blood cells
Immunogen Region	
Specificity	
Immunogen Sequence	