

Anti-Erythroid Beta-Spectrin antibody [DB2] (STJ16100821)

STJ16100821

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Erythroid Beta-Spectrin is suitable for use in Immunohistochemistry, Western Blot and Immunocytochemistry research applications.
Applications	IHC-F/WB/ICC
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	DB2
Concentration	100 ug/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 1% BSA & 0.1% sodium azide
Isotype	IgG1
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Ghost proteins of human red blood cells
Immunogen Region	
Specificity	
Immunogen Sequence	