

Anti-HCAM antibody [DF 1485] (STJ16100826)

STJ16100826

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-HCAM is suitable for use in Immunohistochemistry, Flow Cytometry and Western Blot research applications.
Applications	IHC-F/IHC-P/FC/WB
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	DF 1485
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	IHC-F IHC-P 1:5-1:10 FC WB
Formulation	Contains 0.09% sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	CD44 antigen (PGP-1) purified from lymphocyte membrane
Immunogen Region	
Specificity	
Immunogen	
Sequence	