

Anti-ACTH antibody [57] (STJ16100830)

STJ16100830

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-ACTH is suitable for use in Immunohistochemistry, ELISA and Immunofluorescence research applications.
Applications	IHC-P/ELISA/IF
Host/Source	Mouse
Reactivity	Human/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	57
Concentration	1.0mg/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.09% sodium azide
Isotype	IgG1
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID	72
Gene Symbol	ACTG2
Uniprot ID	ACTH_HUMAN
Immunogen	Synthetic peptide corresponding to the amino acids 1-24 (SYSMEHFRWGKPVGKKRRPVKVYP) of the human ACTH conjugated to Keyhole Limpet Haemocyanin.
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	