

Anti-CD88 antibody [S5/1] (STJ16100838)

STJ16100838

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD88 is suitable for use in Western Blot, Immunohistochemistry and Flow Cytometry research applications.
Applications	WB/IHC-P/FC
Host/Source	Mouse
Reactivity	Human/Rabbit/Ferret/Bovine/Mustelid/Mink

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	S5/1
Concentration	1.0 mg/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.09% sodium azide
Isotype	IgG2a
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	C5aR-peptide: Met1-Asn31.
Immunogen Region	
Specificity	
Immunogen Sequence	



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081