

Anti-Mast Cell Tryptase antibody [AA1] (STJ16100841)

STJ16100841

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Mast Cell Tryptase is suitable for use in Western Blot, Immunohistochemistry, ELISA and Immunofluorescence research applications.
Applications	WB/IHC-P/IHC-F/ELISA/IF
Host/Source	Mouse
Reactivity	Human/Cat/Dog/Monkey/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	AA1
Concentration	1.0 mg/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.09% sodium azide
Isotype	IgG1
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Human mast cell tryptase purified from human lung tissue.
Immunogen Region	
Specificity	
Immunogen Sequence	