

## Anti-CD40 antibody [LOB7/6] (STJ16100844)

STJ16100844

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Mouse monoclonal antibody anti-CD40 is suitable for use in Immunohistochemistry, Flow Cytometry, Immunohistochemistry and Immunoprecipitation research applications.
<b>Applications</b>	IHC-P/FC/IHC-F/IP
<b>Host/Source</b>	Mouse
<b>Reactivity</b>	Human/Dog

### PRODUCT PROPERTIES

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	LOB7/6
<b>Concentration</b>	1.0 mg/mL
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	
<b>Dilution Range</b>	
<b>Formulation</b>	TBS
<b>Isotype</b>	IgG2a
<b>Storage</b>	Store at 2-8°C upon receipt.
<b>Instruction</b>	

### TARGET INFORMATION

<b>Gene ID</b>	958
<b>Gene Symbol</b>	CD40
<b>Uniprot ID</b>	TNR5_HUMAN
<b>Immunogen</b>	CD40 Ig (Fc) fusion protein containing the EC region of human CD40 and Fc region of human IgG.
<b>Immunogen Region</b>	
<b>Specificity</b>	
<b>Immunogen Sequence</b>	

