

Anti-MyeloHRP antibody [2C7] (STJ16100847)

STJ16100847

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-MyeloHRP is suitable for use in Western Blot, Immunohistochemistry, Flow Cytometry, Immunohistochemistry and ELISA research applications.
Applications	WB/IHC-P/FC/IHC-F/ELISA
Host/Source	Mouse
Reactivity	Human/Dog

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	2C7
Concentration	1.0 mg/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.09% sodium azide
Isotype	IgG1
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Human myeloperoxidase
Immunogen Region	
Specificity	
Immunogen Sequence	