

Anti-CD21 antibody [2G9] (STJ16100848)

STJ16100848

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD21 is suitable for use in Immunohistochemistry and Immunohistochemistry research applications.
Applications	IHC-P/IHC-F
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	2G9
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.09% sodium azide
Isotype	IgG2a
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Recombinant human CD21.
Immunogen Region	
Specificity	
Immunogen	
Sequence	