

Anti-CD68 antibody [514H12] (STJ16100849)

STJ16100849

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD68 is suitable for use in Immunohistochemistry and Immunofluorescence research applications.
Applications	IHC-P/IHC-F/IF
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	514H12
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.09% sodium azide
Isotype	IgG2a
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	968
Gene Symbol	CD68
Uniprot ID	CD68_HUMAN
Immunogen	Fusion protein corresponding to external domain of human CD68.
Immunogen Region	
Specificity	
Immunogen	
Sequence	