

## Anti-Mcm2 antibody [CRCT2.1] (STJ16100851)

STJ16100851

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Mouse monoclonal antibody anti-Mcm2 is suitable for use in Western Blot and Immunohistochemistry research applications.
<b>Applications</b>	WB/IHC-P
<b>Host/Source</b>	Mouse
<b>Reactivity</b>	Human

### PRODUCT PROPERTIES

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	CRCT2.1
<b>Concentration</b>	1.0mg/mL
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	
<b>Dilution Range</b>	
<b>Formulation</b>	PBS with 0.09% sodium azide
<b>Isotype</b>	IgG1
<b>Storage Instruction</b>	Store at 2-8°C upon receipt.

### TARGET INFORMATION

<b>Gene ID</b>	<a href="#">4171</a>
<b>Gene Symbol</b>	<a href="#">MCM2</a>
<b>Uniprot ID</b>	<a href="#">MCM2_HUMAN</a>
<b>Immunogen</b>	Bacterially expressed human Mcm2.
<b>Immunogen Region</b>	
<b>Specificity</b>	
<b>Immunogen Sequence</b>	