

Anti-Actin Alpha-muscle antibody [HHF35] (STJ16100854)

STJ16100854

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Actin Alpha-muscle is suitable for use in EM, ELISA, Immunocytochemistry, Immunohistochemistry and Western Blot research applications.
Applications	EM/ELISA/ICC/IHC-F/IHC-P/WB
Host/Source	Mouse
Reactivity	Human/Chicken/Monkey/Rabbit/Rat/Swine/Zebrafish

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	HHF35
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.09% sodium azide
Isotype	IgG1
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	SDS-extracted protein fraction from Human myocardium
Immunogen Region	
Specificity	
Immunogen Sequence	



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081