

Anti-Aquain 1 antibody [1/A5F6] (STJ16100858)

STJ16100858

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Aquain 1 is suitable for use in Western Blot, Immunohistochemistry, ELISA and Immunofluorescence research applications.
Applications	WB/IHC-P/IHC-F/ELISA/IF
Host/Source	Mouse
Reactivity	Human/Rabbit/Rat/Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	1/A5F6
Concentration	1.0 mg/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.1% sodium azide and 0.1% BSA
Isotype	IgG1
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Synthetic peptide corresponding to amino acids 249-269 of aquaporin 1.
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	