

Anti-Granzyme A antibody [GA6] (STJ16100861)

STJ16100861

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Granzyme A is suitable for use in Western Blot and Immunohistochemistry research applications.
Applications	WB/IHC-P
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	GA6
Concentration	1.0 mg/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.09% sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	3001
Gene Symbol	GZMA
Uniprot ID	GRAA_HUMAN
Immunogen	Recombinant human Granzyme A.
Immunogen Region	
Specificity	
Immunogen	
Sequence	