

Anti-TGF Beta antibody [TB21] (STJ16100866)

STJ16100866

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-TGF Beta is suitable for use in Western Blot, Immunohistochemistry, ELISA, Flow Cytometry and Immunohistochemistry research applications.
Applications	WB/IHC-P/ELISA/FC/IHC-F
Host/Source	Mouse
Reactivity	Human/Mink/Sheep/Rabbit/Rat/Mustelid/Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	TB21
Concentration	1.0 mg/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS
Isotype	IgG1
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Human Transforming Growth Factor Beta 1 from platelets.
Immunogen Region	
Specificity	
Immunogen Sequence	