

**Anti-Macrophages antibody [MAC387] (STJ16100868)**

STJ16100868

**GENERAL INFORMATION**

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Mouse monoclonal antibody anti-Macrophages is suitable for use in Immunohistochemistry, Flow Cytometry and Immunohistochemistry research applications.
<b>Applications</b>	IHC-P/FC/IHC-F
<b>Host/Source</b>	Mouse
<b>Reactivity</b>	Human/Rat/Rabbit/Pygmy Hippopotamus/Pig/Rhesus/Mink/Marmoset/Bovine/Baboon/Cat/Cynomolgus Monkey/Fallow Deer/Goat/Dog/Guinea Pig/Horse

**PRODUCT PROPERTIES**

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	MAC387
<b>Concentration</b>	1.0 mg/mL
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	
<b>Dilution Range</b>	
<b>Formulation</b>	PBS with 0.09% sodium azide
<b>Isotype</b>	IgG1
<b>Storage</b>	Store at 2-8°C upon receipt.
<b>Instruction</b>	

**TARGET INFORMATION**

<b>Gene ID</b>	
<b>Gene Symbol</b>	
<b>Uniprot ID</b>	
<b>Immunogen</b>	Human monocytes.
<b>Immunogen</b>	
<b>Region</b>	
<b>Specificity</b>	
<b>Immunogen</b>	
<b>Sequence</b>	