

**Anti-SKP2 antibody [SKP2-8D9] (STJ16100869)**

STJ16100869

**GENERAL INFORMATION**

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Mouse monoclonal antibody anti-SKP2 is suitable for use in Western Blot, Immunohistochemistry, Immunocytochemistry, Immunofluorescence, ELISA and Immunoprecipitation research applications.
<b>Applications</b>	WB/IHC/ICC/IF/ELISA/IP
<b>Host/Source</b>	Mouse
<b>Reactivity</b>	Human/Non-Human Primates

**PRODUCT PROPERTIES**

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	SKP2-8D9
<b>Concentration</b>	0.5 mg/mL
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	
<b>Dilution Range</b>	WB IHC ICC IF 1:250 ELISA IP
<b>Formulation</b>	PBS with 0.1% sodium azide
<b>Isotype</b>	IgG1k
<b>Storage Instruction</b>	Store at 2-8°C upon receipt.

**TARGET INFORMATION**

<b>Gene ID</b>	6502
<b>Gene Symbol</b>	SKP2
<b>Uniprot ID</b>	SKP2_HUMAN
<b>Immunogen</b>	Recombinant, full-length human p45SKP2 protein
<b>Immunogen Region</b>	
<b>Specificity</b>	
<b>Immunogen Sequence</b>	