

Anti-CD63 antibody [MEM-259] {PE} (STJ16100872)

STJ16100872

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD63 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	MEM-259
Concentration	
Conjugation	PE
Purification	
Dilution Range	
Formulation	Stabilizing phosphate buffered saline (PBS) solution containing 15 mM sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	967
Gene Symbol	CD63
Uniprot ID	CD63_HUMAN
Immunogen	HPB-ALL T cell line
Immunogen Region	
Specificity	
Immunogen Sequence	