

Anti-NCAM antibody [123C3] (STJ16100875)

STJ16100875

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD56 is suitable for use in Immunohistochemistry and Immunofluorescence research applications.
Applications	IHC-P/IHC-F/IF
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	123C3
Concentration	100 ug/mL
Conjugation	Unconjugated
Purification	
Dilution Range	IF 1:50 IHC 1:50
Formulation	PBS with 0.1% BSA and 0.02% sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen Region
Specificity
Immunogen
Sequence



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081