

Anti-Carcinoembryonic Antigen antibody [PARLAM4] (STJ16100878)

STJ16100878

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Carcinoembryonic antigen is suitable for use in Immunohistochemistry and Western Blot research applications.
Applications	IHC-F/IHC-P/WB
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	PARLAM4
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	IHC-F IHC-P 1:5-1:10 WB
Formulation	Contains 0.09% sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	1048
Gene Symbol	CEACAM5
Uniprot ID	CEAM5_HUMAN
Immunogen	Isolated from human colon carcinoma cells
Immunogen Region	
Specificity	
Immunogen Sequence	