

Anti-NCAM antibody [NKI-nbl-1] (STJ16100882)

STJ16100882

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD56 is suitable for use in Immunohistochemistry and Immunohistochemistry research applications.
Applications	IHC-P/IHC-F
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	NKI-nbl-1
Concentration	100 ug/mL
Conjugation	Unconjugated
Purification	
Dilution Range	IHC-P IHC-F 1:10
Formulation	PBS with 0.1% BSA and 0.02% sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen Region
Specificity
Immunogen Sequence