

Anti-Synaptophysin antibody [SY38] (STJ16100883)

STJ16100883

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Synaptophysin is suitable for use in Immunohistochemistry, Immunofluorescence and Western Blot research applications.
Applications	IHC-F/IHC-P/IF/WB
Host/Source	Mouse
Reactivity	Bovine/Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	SY38
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	IHC-F IHC-P 1:10 IF WB
Formulation	PBS with 0.5% BSA and 0.09% sodium azide
Isotype	IgG1
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID	6855
Gene Symbol	SYP
Uniprot ID	SYPH_HUMAN
Immunogen	Synaptophysin from presynaptic vesicles, prepared from bovine brain
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	