

**Anti-P110 glycoprotein antibody [LRP-56] (STJ16100886)**

STJ16100886

**GENERAL INFORMATION**

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Mouse monoclonal antibody anti-P110 glycoprotein is suitable for use in Flow Cytometry, Immunocytochemistry, Immunohistochemistry and Immunohistochemistry research applications.
<b>Applications</b>	FC/ICC/IHC-F/IHC-P
<b>Host/Source</b>	Mouse
<b>Reactivity</b>	Human

**PRODUCT PROPERTIES**

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	LRP-56
<b>Concentration</b>	100 ug/mL
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	
<b>Dilution Range</b>	FC 1:50 ICC IHC-F 1:20-1:50 IHC-P 1:20
<b>Formulation</b>	supernatant with 1% BSA and 0.1% sodium azide
<b>Isotype</b>	IgG2b
<b>Storage Instruction</b>	Store at 2-8°C upon receipt.

**TARGET INFORMATION**

<b>Gene ID</b>
<b>Gene Symbol</b>
<b>Uniprot ID</b>
<b>Immunogen</b>
<b>Immunogen Region</b>
<b>Specificity</b>
<b>Immunogen Sequence</b>