

**Anti-LRP/MVP antibody [LMR5] (STJ16100893)**

STJ16100893

**GENERAL INFORMATION**

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Mouse monoclonal antibody anti-LRP/MVP is suitable for use in Immunohistochemistry and Flow Cytometry research applications.
<b>Applications</b>	IHC-F/IHC-P/FC
<b>Host/Source</b>	Mouse
<b>Reactivity</b>	Human

**PRODUCT PROPERTIES**

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	LMR5
<b>Concentration</b>	250 ug/mL
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	
<b>Dilution Range</b>	
<b>Formulation</b>	supernatant with 1% BSA and 0.1% sodium azide
<b>Isotype</b>	IgG
<b>Storage Instruction</b>	Store at 2-8°C upon receipt.

**TARGET INFORMATION**

<b>Gene ID</b>	9961
<b>Gene Symbol</b>	MVP
<b>Uniprot ID</b>	MVP_HUMAN
<b>Immunogen</b>	
<b>Immunogen Region</b>	
<b>Specificity</b>	
<b>Immunogen</b>	
<b>Sequence</b>	