

Anti-MRP5 antibody [M5I-1] (STJ16100899)

STJ16100899

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-MRP5 is suitable for use in Immunocytochemistry, Immunohistochemistry, Flow Cytometry and Western Blot research applications.
Applications	ICC/IHC-F/FC/WB
Host/Source	Rat
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	M5I-1
Concentration	250 ug/mL
Conjugation	Unconjugated
Purification	
Dilution Range	ICC 1:20-1:50 IHC-F 1:20 FC WB 1:50
Formulation	supernatant with 0.7% BSA and 0.1% sodium azide
Isotype	IgG2a
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	10057
Gene Symbol	ABCC5
Uniprot ID	MRP5_HUMAN
Immunogen	MRP5 aa 82-168
Immunogen Region	
Specificity	
Immunogen Sequence	