

Anti-LRP/MVP antibody [MVP-37] (STJ16100902)

STJ16100902

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-LRP/MVP is suitable for use in Immunohistochemistry, Western Blot and Flow Cytometry research applications.
Applications	IHC-F/IHC-P/WB/FC
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	MVP-37
Concentration	250 ug/mL
Conjugation	Unconjugated
Purification	
Dilution Range	IHC-F IHC-P 1:20-1:50 WB 1:50 FC
Formulation	supernatant with 0.7% BSA and 0.1% sodium azide
Isotype	IgG2b
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	9961
Gene Symbol	MVP
Uniprot ID	MVP_HUMAN
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	