

Anti-Breast Cancer Resistance Protein antibody [BXP-34] (STJ16100903)

STJ16100903

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Breast Cancer Resistance Protein is suitable for use in Immunohistochemistry and Flow Cytometry research applications.
Applications	IHC-F/FC
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	BXP-34
Concentration	250 ug/mL
Conjugation	Unconjugated
Purification	
Dilution Range	IHC-F 1:10-1:20 FC
Formulation	PBS with 1% BSA and 0.1% sodium azide
Isotype	IgG1
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID	9429
Gene Symbol	ABCG2
Uniprot ID	ABCG2_HUMAN
Immunogen	Mitoxanthrone resistant, BCRP overexpressing cell line MCF7 MR
Immunogen Region	
Specificity	
Immunogen Sequence	