

Anti-P193/VPARP antibody [P193-4] (STJ16100905)

STJ16100905

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-P193/VPARP is suitable for use in Immunohistochemistry, Western Blot and Flow Cytometry research applications.
Applications	IHC-P/IHC-F/WB/FC
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	P193-4
Concentration	250 ug/mL
Conjugation	Unconjugated
Purification	
Dilution Range	IHC-P 1:10-1:20 IHC-F WB FC
Formulation	PBS with 1% BSA and 0.1% sodium azide
Isotype	IgG2b k
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	VPARP aa 408-611
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	