

**Anti-P193/VPARP antibody [P193-6] (STJ16100906)**

STJ16100906

**GENERAL INFORMATION**

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Mouse monoclonal antibody anti-P193/VPARP is suitable for use in Immunohistochemistry, Western Blot and Flow Cytometry research applications.
<b>Applications</b>	IHC-P/IHC-F/WB/FC
<b>Host/Source</b>	Mouse
<b>Reactivity</b>	Human

**PRODUCT PROPERTIES**

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	P193-6
<b>Concentration</b>	250 ug/mL
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	
<b>Dilution Range</b>	IHC-P 1:10-1:20 IHC-F WB FC
<b>Formulation</b>	PBS with 1% BSA and 0.1% sodium azide
<b>Isotype</b>	IgG2b k
<b>Storage</b>	Store at 2-8°C upon receipt.
<b>Instruction</b>	

**TARGET INFORMATION**

<b>Gene ID</b>	
<b>Gene Symbol</b>	
<b>Uniprot ID</b>	
<b>Immunogen</b>	VPARP aa 408-611
<b>Immunogen</b>	
<b>Region</b>	
<b>Specificity</b>	
<b>Immunogen</b>	
<b>Sequence</b>	