

Anti-MVP antibody [1014] (STJ16100910)

STJ16100910

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Major Vault Protein is suitable for use in Western Blot, Immunofluorescence, Immunohistochemistry and Immunoprecipitation research applications.
Applications	WB/IF/IHC/IP
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	1014
Concentration	Can be provided as 100 µg/mL, 500 µg/mL or 1mg/mL.
Conjugation	Unconjugated
Purification	Affinity purified from tissue culture.
Dilution Range	Immunoblot (1-2 µg/ml). Immunofluorescence (0, 5-1, 0 µg/ml). Immunohistology (1-2 µg/ml for 30-60 minutes at RT; information on a suitable antigen retrieval method for staining of formalin-fixed tissues is unavailable to date). Immunoprecipitation (
Formulation	PBS with 0.02% Sodium Azide.
Isotype	IgG1k
Storage	Store for up to 1 year at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID	9961
Gene Symbol	MVP
Uniprot ID	MVP_HUMAN
Immunogen	A BALB/c mouse was immunized with affinity purified nuclear extract proteins. Fusion partner: Sp2/O-Ag14.
Immunogen Region	
Specificity	
Immunogen Sequence	

