

Anti-MRP6 antibody [M6II-21] (STJ16100914)

STJ16100914

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-MRP6 is suitable for use in Immunocytochemistry, Immunohistochemistry, Flow Cytometry and Western Blot research applications.
Applications	ICC/IHC-F/FC/WB
Host/Source	Rat
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	M6II-21
Concentration	100 ug/mL
Conjugation	Unconjugated
Purification	
Dilution Range	ICC IHC-F 1:20-1:50 FC WB 1:20-1:50
Formulation	supernatant with 1% BSA and 0.1% sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	368
Gene Symbol	ABCC6
Uniprot ID	MRP6_HUMAN
Immunogen	ABCC6 aa 764-964
Immunogen Region	
Immunogen Specificity	
Immunogen Sequence	