

**Anti-MRP6 antibody [M6II-31] (STJ16100915)**

STJ16100915

**GENERAL INFORMATION**

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Rat monoclonal antibody anti-MRP6 is suitable for use in Immunocytochemistry, Immunohistochemistry, Flow Cytometry and Western Blot research applications.
<b>Applications</b>	ICC/IHC-F/IHC-P/FC/WB
<b>Host/Source</b>	Rat
<b>Reactivity</b>	Human

**PRODUCT PROPERTIES**

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	M6II-31
<b>Concentration</b>	50 ug/mL
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	
<b>Dilution Range</b>	ICC IHC-F 1:20-1:50 IHC-P FC WB 1:20-1:50
<b>Formulation</b>	supernatant with 0.7% BSA and 0.1% sodium azide
<b>Isotype</b>	IgG2a
<b>Storage</b>	Store at 2-8°C upon receipt.
<b>Instruction</b>	

**TARGET INFORMATION**

<b>Gene ID</b>	<a href="#">368</a>
<b>Gene Symbol</b>	<a href="#">ABCC6</a>
<b>Uniprot ID</b>	<a href="#">MRP6_HUMAN</a>
<b>Immunogen</b>	ABCC6 aa 764-964
<b>Immunogen Region</b>	
<b>Specificity</b>	
<b>Immunogen Sequence</b>	