

Anti-Breast Cancer Resistance Protein antibody [BXP-9] (STJ16100916)

STJ16100916

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-Breast Cancer Resistance Protein is suitable for use in Immunocytochemistry, Immunohistochemistry, Western Blot and Flow Cytometry research applications.
Applications	ICC/IHC-F/WB/FC
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	BXP-9
Concentration	250 ug/mL
Conjugation	Unconjugated
Purification	
Dilution Range	ICC 1:20-1:50 IHC-F 1:20 WB 1:50 FC
Formulation	cell culture supernatant with 0.7% BSA and 0.1% sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	9429
Gene Symbol	ABCG2
Uniprot ID	ABCG2_HUMAN
Immunogen	BCRP aa 221-394
Immunogen Region	
Immunogen Specificity	
Immunogen Sequence	