

Anti-CD13 antibody [WM15] (STJ16100918)

STJ16100918

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD13 is suitable for use in Flow Cytometry, Immunoprecipitation and Immunohistochemistry research applications.
Applications	FC/IP/IHC-F
Host/Source	Mouse
Reactivity	Human/Non-Human Primates

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	WM15
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	Phosphate buffered saline (PBS) solution with 15 mM sodium azide
Isotype	IgG1
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID	290
Gene Symbol	ANPEP
Uniprot ID	AMPN_HUMAN
Immunogen	Human AML cells
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	