

Anti-MRP4 antibody [M4I-10] (STJ16100921)

STJ16100921

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-MRP4 is suitable for use in Immunocytochemistry, Immunohistochemistry, Western Blot and Flow Cytometry research applications.
Applications	ICC/IHC-F/WB/FC
Host/Source	Rat
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	M4I-10
Concentration	100 ug/mL
Conjugation	Unconjugated
Purification	Protein G purified
Dilution Range	ICC 1:20-1:50 IHC 1:20 WB 1:20-1:50
Formulation	Tissue culture supernatant with 0.1% BSA and 0.09% sodium azide
Isotype	IgG2a
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID	10257
Gene Symbol	ABCC4
Uniprot ID	MRP4_HUMAN
Immunogen	Maltose Binding Protein and aa 372-431 of MRP4
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	